



HCP-7000 CHART PROJECTOR

- Semi-permanent lifetime LED light source
- Brighter and clearer than conventional bulb light
- Attractive Streamlined Shape and Blue Indicator
- RC-R Type (Russian Chart) added



SPECIFICATIONS

| | |
|--------------------------|--|
| Chart | 41 charts 34 masks, Red/Green & Polarization Filters |
| Projection Distance | 2.5~8m |
| Chart Rotation Speed | Average 0.15 sec |
| Projection Magnification | 30x at 5m |
| Power Saving | Automatic Switch off (10 min) |
| Program | 2 programs with a maximum of 30 charts each |
| Tilt Angle | 15° |
| Power Supply | 100-120V 50Hz : 0.6A, 200-240V 60Hz : 0.3A |
| Lamp | LED 4W |
| Dimension | 270(W) x 182(D) x 230(H)mm |
| Weight | 3.44kg |

Designs and details can be changed without prior notice for the purposes of improvement.

Huvitz

HUVITZ Co., Ltd. 38, Burim-ro 170beon-gil, Dongan-gu,
Anyang-si, Gyeonggi-do, 14055, Republic of Korea
Tel:+82-31-442-8868 Fax:+82-31-477-8617 <http://www.huvitz.com>

E2ARAA-16-00001-2-170426



HCP-7000 CHART PROJECTOR

Huvitz Re:define. Re+create