OCULUS | BIOM® ready

Single-use wide angle viewing system






OCULUS BIOM® ready

The world's first single-use wide angle viewing system.

From the inventors of non-contact wide angle viewing, comes the worlds first single-use wide angle viewing system. The BIOM® ready provides the retinal surgeon the same unparalled view they have come to rely on for over 25 years, while addressing the need for maximizing operating room efficiency.

- Perfect balance between flexibility and high optical quality
- Uncompromised performance you have come to expect from the OCULUS BIOM®
- No sterilization "down time" for increased O.R. efficiency



Incorporates the new BIOM® HD Disposable Lens for unparalleled visual clarity


- Excellent depth of field for better stereopsis
- HD clarity under high magnification reduces the need for a contact lens
- Full-field clarity for decreased scleral depressing during panretinal laser
- Improved view during air/fluid exchanges



130° wide angle field of view with outstanding resolution in the periphery



New reduction lens is designed to optimize the view of the retina and the limbus, increasing surgical efficiency



HD clarity under high magnification reduces the need for a contact lens





Always ready, always a perfect view.







Connects easily to your microscope



Ready for use

Specifically designed to optimize performance with either a f=175 mm or f=200 mm microscope objective lens.

54413: BIOM® ready, optimized for f = 175 mm 54414: BIOM® ready, optimized for f = 200 mm

Adaptations

SDI® 4/BIOM® ready adaptations are currently available for the following microscope types:	
Alcon	LuxOR
Kaps	SOM
Leica	M500/M501, M620, M650, M690, M841, M820, M822, M844
Möller	Ophtamic 900, Hi-R 900, EOS 900, Allegra 900
Takagi	OM-18, OM-8
Topcon	OMS-110, OMS-600, OMS-610, OMS-650, OMS-710, OMS-800 Standard, OMS-800 Pro, OMS-850 Standard, OMS-850 Pro
Zeiss	OPMI VISU 150/160, VISU 200/210, Lumera 700, Lumera T, Lumera i, OPMI CS with Retrolux 1/3/CS, OPMI MDI/MDO/MDU, OPMI 1/6

(← in accordance with Medical Device Directive 93/42/EEC

WWW.OCULUSSURGICAL.COM

Headquarters:

OCULUS Optikgeräte GmbH Postfach • 35549 Wetzlar • GERMANY Tel. +49-641-2005-0 • Fax +49-641-2005-295 Email: export@oculus.de • www.oculus.de

OCULUS Surgical, Inc.

562 NW Mercantile Place, Ste. 104 Port St. Lucie • FL 34986 • USA Toll free 855 SDI BIOM • Tel. 772 236 2622 Email: info@oculussurgical.com • www.oculussurgical.com

