





COSTRUZIONE
STRUMENTI
OFTALMICI



MIZAR
VIDEOCAMERA DIGITALE
5 MPIXEL

MIZAR
DIGITAL CAMERA
5 MPIXEL



per Lampada a Fessura SL9800 | SL9900 | SL9900 Elite
for Slit Lamps SL9800 | SL9900 | SL9900 Elite

VIDEOCAMERA DIGITALE MIZAR 5 MPIXEL

MIZAR può essere applicata alle lampade a fessura SL9800, SL9900 e SL9900 Elite con microscopio 3X, 5X e Zoom. Questa videocamera si distingue per le sue elevate prestazioni ed è caratterizzata da un'ottima resa cromatica. Per mezzo del software applicativo Phoenix, durante l'esame biomicroscopico, consente la visualizzazione dei più piccoli dettagli, l'acquisizione di immagini e filmati ad alta risoluzione e la loro l'archiviazione.

DIGITAL CAMERA MIZAR 5 MPIXEL

MIZAR can be applied to SL9800, SL9900 and SL9900 Elite slit lamps, with 3X, 5X and Zoom microscope. This video camera stands out for its high performance and is characterized by excellent color rendering. While the biomicroscopic examination, the Phoenix application software, allows the display of the smallest details and the acquisition of high-resolution images and videos, their storage and management for further review.



CARATTERISTICHE TECNICHE		TECHNICAL DATA	
Risoluzione	5Mpx	Resolution	5Mpx
Risoluzione pixel	2448x2048	Pixel resolution	2448x2048
Connessione	USB 3.0	Camera interface	USB 3.0
Frame Rate	35 fps	Frame Rate	35 fps

CPU: I3 o superiore (consigliata I5), RAM: 4GB (consigliati 8GB), OS: Win10x64, Win8x64, Win7x32, Win7x64, GPU: 1GB memoria condivisa, Risoluzione schermo: 1280x960 o superiore

CPU: I3 or higher (suggested I5), RAM: 4 GB or higher (suggested 8 GB), OS: Win10x64, Win8x64, Win7x32, Win7x64, GPU: 1 Gbyte not shared, RESOLUTION: 1280x960 or higher



COSTRUZIONE STRUMENTI OFTALMICI
Via degli Stagnacci 12/E
50018 - Scandicci (FI) | Italy
tel: +39 055 722191 | email cs@csitalia.it
WWW.CSOITALIA.IT | WWW.CSOEVENTI.IT

CO104 | Rev. 10/03/2022