



Portable Slit Lamp

SL1911D

High Quality Optical Design Inside

Listen to the clinical

Features:

- Multiple coating glass deliver high quality optics
- Warm white **LED** similar to halogen bulb with **20000** hours life span
- Continuous working time up to **6** hours
- **1mm square light patch** for assessing anterior chamber flare
- Two magnification of 10X and 16X switched by lever

Specification

MICROSCOPE

Optics	Converging binoculars at 13°
Magnification	10X and 16X
Field of View	φ10mm at 16X magnification φ15mm at 10X magnification
PD Range	50 ~ 74 mm
Diopter Range	-7D ~ +7D

ILLUMINATION

Light Source	Warm white LED
Slit Length	12mm
Slit Width	0.15mm, 0.5mm, 0.8mm, 1.6mm, 12mm circle, 1mm square
Filter	Red free, Cobalt blue, Neutral density
Slit Angle	+/- 60°

SL1911D

POWER

Battery	7.4V Rechargeable Li-ion battery, 1600mAh Working time: About 6 hours (at maximum light intensity)
Charging Station	Input voltage: AC 100-240V, 50-60Hz Out voltage: DC 12V, 2A

OTHERS

Dimension	124(W) x 195(D) x 214(H)mm
Weight	about 894g (with battery)
Accessories	Aluminum carrying case, Battery(1pc), Charging station, Power adaptor, Forehead support, Diopter adjusting bar, User manual

 FDA Approval

Shanghai VisuScience Meditech Co.,Ltd.

Add: No. 344 Sanlin Road, Pudong New Area, Shanghai, China

Tel: +86-21- 34973659 | Fax: +86-21- 34973659

E-mail: info@visuscience.com

Follow us on



www.visuscience.com